



An ISO 9001:2015 Company
An ISO 13485:2012 Company

Scopy Camera (EC01)



MEDICAM introduces the new state of art scopy camera which can be used with almost all the standard fibre optic scopes.

Features rich camera with features like manual white balance on the camera head, ergonomic design, and lightweight body makes the camera user friendly in true sense.

Features

- Ergonomic design
- Lightweight
- High Resolution
- Manual white balance on the Camera Head
- Light weight control box
- Double A/V out
- Fits on any coupler
- Easy to handle

Specifications

- Sensor – 1/3 Sony CCD
- Horizontal Resolution – 650TVL
- Total Pixels – 1020 (H) X 596 (V) 610 K
- Min. illumination – 0.1 Lux
 - Output – 1vp-p composite [OHMS]
 - Output Ports – 02 nos
 - Power – 12V DC
 - Shutter Speed- 1/50 ~ 90,000(PAL)

Dimensions

- Camera Head – Ø 31.50mm
- Camera Length – 62.80mm
- Control Box – 165 (L) X 130 (W) X 62 (H)

Operating Condition

- Temperature Recommendation - 5°C ~ 40°C



Sh. No. 20/22, Bhivandiwala Bldg. Dr. Ambedkar Road, Kalachowki, Mumbai - 400033
Email: info@hemantmedicam.com, sales@medicam.co.in, hemantmedicam@yahoo.co.in
Ph No.: +91 22 24700872, 8767303439
Website: www.medicam.in | www.hemantmedicam.com | www.medicam.co.in